Leadex 70 AC

Intraoraler Röntgenstrahler - intraoral x-ray system

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OPTIONAL

- ▶ SECOND CONTROL BUTTON
- ▶ EXTERNAL SIGNAL LAMP
- ▶ EXTERNAL CONTROL BUTTON

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Fig. 2A, 2B, 2C give the overall dimensions of the possible supply conditions:

- ▶ **WANDARM 400:** Länge 41cm – 16,2"
- ▶ **WANDARM 800:** Länge 82,5cm – 32,5"
- ▶ **WANDARM 1100:** Länge 110cm – 43,5"

- ▶ **BRACKET 400:** length 41cm – 16,2"
- ▶ **BRACKET 800:** length 82,5cm – 32,5"
- ▶ **BRACKET 1100:** length 110cm – 43,5"

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Fig. 3 and 4 show the typical dimensions of the radiographic system:

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WALL WORKS PREPARATION

- 1) Choose on the wall the fixing position referring the working axis and taking into consideration the electric system type (enclosed, outside or flying type) according to the position assigned to the control board.
- 2) Mark on the wall the holes position using the carton template and check its perpendicularity with a plumb line. If the electric system is enclosed, mark also the corresponding hole.
- 3) Make 6 holes on the wall starting with a $\varnothing 7$ drill and gradually enlarging them up to $\varnothing 12$. This procedure is made to avoid the demolition of the white coat and to keep the distance between centers under control. When the walls are made of hollow, full or cement bricks use $\varnothing 12$ metallic dowels provided with grub screw which cannot be lost and with M 6 separate screw with hexagonal head and washer. With not reliable walls it is necessary to make a suitable reinforcement.
- 4) Mount the set approaching it parallelly to the wall and tighten the six screws alternatively; if the wall is not perfectly plane interpose some washers.
- 5) When the assembly and the dynamic testing of the equipment is made, check that the screws are perfectly tightened to eliminate any clearance between the wall and the plate.

NON USARE TASSELLI DI MATERIALE PLASTICO
DO NOT USE PLASTIC DOWELS
NE PAS UTILISER DE TAQUETS FABRIQUES EN MATERIEL PLASTIQUE

345 mm

COD. 0.06.00435

minimum

cm 130

110 mm

Major Dental Equipment Utility Specifications

Ritter Leadex 70 DC Intra Oral X Ray

Electrical

- | | | | |
|----|---|----------------------|---|
| 1. | P | POWER | <p>supply GPO <i>minimum</i> 240v/10a power to marked location area (refer x- ray mounting template), Leave 500mm tail as requested.</p> <p>Power Supply must be installed in accordance with AS/NZS 3003.2011 for Body Protection with R.C.D (Residual Current Device)</p> <p><i>Check x ray unit for Electrical circuit requirements.</i></p> <p><i>Should be wired to same REMOTE isolation on/off switch as dental chair system.</i></p> |
| 2 | R | REMOTE | <p>for remote control connection, supply 3 core SWITCHING flex cable from x-ray location to selected area on wall. leave tail of 600mm each end.</p> |
| 3 | S | STRUCTURAL REINFORCE | <p>Where required, ensure wall mounting point for x-ray is reinforced with minimum 700x700x32mm board (spanning/mounted to internal wall studs) at location for mounting of x-ray unit.</p> <p>Bottom of board to be 110cm above floor level</p> <p>Reinforcement of wall should be made internally to wall prior to Gyprock sheeting – ensure no gap between reinforcement board and Gyprock sheeting.</p> <p>Refer Manufacturer's template for exact location for reinforcement of wall for x-ray in relation to floor and dental system.</p> |