

Breathe through your nose.

ClearView

SINGLE-USE NASAL HOODS



Now you can actually see your patient breathing through their nose.

Making a Difference

Designed for the comfort and safety of your patient and enhancement of your clinical treatment, Accutron's ClearView redefines the nitrous oxide nasal hood market.

- Clear outer hood provides a visual of patient respiration – condensation appears when a patient breathes through their nose.
- Low profile expands clinical field of view and allows for easy access to the mouth.
- Soft inner hood creates desirable nasal/facial seal reducing seepage of exhaled gas.
- Hood-in-hood design – outer hood provides additional assurance that exhaled gases are collected minimizing ambient nitrous oxide in the dental operator.
- Colorful and lightly scented inner hood engages and delights patients of all ages.
- Latex-free material protects patients and staff who have latex allergies.
- Single-use and individually wrapped nasal hoods offer safety assurance; saves staff time with operator turnaround.
- Pedito and adult sizes available in a variety of scents plus unscented.



Latex-free

PEDO SIZES:

- Birthday Bubblegum
- Outlaw Orange
- Groovy Grape
- Unscented
- Sassy Strawberry

ADULT SIZES:

- Unscented
- Outlaw Orange
- French Vanilla
- Fresh Mint
- Birthday Bubblegum



Scan the QR code with a QR Reader from your smart phone or mobile device for a direct link to the **free sample** request page.

**For FREE samples call 800-531-2221
or go to www.accutron-inc.com**

*ClearView is available through your
authorized dental distributor.*



Relentless Innovation. Creative Solutions.

ClearView

SINGLE-USE NASAL HOODS

Safety First

Accutron scavenging circuits help to secure a safe and productive work environment for patients, doctor and staff.

- X-ray compatible circuitry saves time if images are required during treatment.
- Easy-to-clean, autoclavable circuitry for protection against cross-contamination.
- Components and rubber goods are latex-free for safety all around.
- Made in the U.S.A. for maximum quality standards.
- Bright white finish complements clinical operatory environments.
- Flexible system accommodates most brands of flowmeters.

Remember to check your scavenging circuit's breathing bag and corrugated tubing for tears and holes on a regular basis.

ClearView Scavenging Circuits are available through your authorized dental distributor.

ClearView Scavenging Circuit II (43003)



ClearView Standard Scavenging Circuit

43003 ClearView Scavenging Circuit II

Includes scavenging circuit, rubber goods (breathing bag and corrugated tube), and a complementary sampling of ClearView Single-Use Nasal Hoods

43005 ClearView Scavenging Circuit I

Includes scavenging circuit and a complementary sampling of ClearView Single-Use Nasal Hoods

ClearView Scavenging Circuit for Remote Flow System (RFS) and Accu-Vac

43004 ClearView Scavenging Circuit II

Includes scavenging circuit, rubber goods (breathing bag and corrugated tube), and a complementary sampling of ClearView Single-Use Nasal Hoods

43006 ClearView Scavenging Circuit I

Includes scavenging circuit and a complementary sampling of ClearView Single-Use Nasal Hoods



1733 Parkside Lane • Phoenix, AZ 85027
Toll-free: 800.531.2221 • Local Phone: 623.780.2020 • Fax: 623.780.0444

www.accutron-inc.com

Part No. 27689 - Rev. 00 - 1/2012

Breathe through your nose. Breathe through your nose. Breathe through your nose.