



SPITTOON VALVE

Spittoon valve for draining the waste from the spittoon bowl off into the suction hose of a central suction system. With integrated prefilter and automatic actuation of the suction engine.

Spittoon valve

Important product characteristics:

- Easy installation in many treatment units due to minimum dimensions and flexible connections.
- Prevents suction noise at the basin through a hermetically sealed closing membrane.
- Easily accessible prefilter prevents bigger parts rinsed into the basin.
- Automatic actuation of the suction system as required.
- For cleaning purposes, the outlet pipe of the spittoon bowl may be manually sucked empty using a probe.

Technical Data:

- | | |
|-----------------------------|---------------------------|
| ➤ Power supply: | 24 V AC or 24 V DC |
| ➤ Frequency: | 50/60 Hz |
| ➤ Max. current consumption: | 100 mA |
| ➤ Max. power input: | 3 VA |
| ➤ Low pressure range: | 80 mbar - 250 mbar |
| ➤ Possible suction system: | wet or dry vacuum systems |
| ➤ Max. water flow rate: | 3,5 l/min |
| ➤ Dimensions (H x W x D): | 125 x 100 x 120 mm |

Scope of delivery:

1 complete spittoon valve

1 manual for installation, operation and maintenance

Spittoon valve

Order no.: 40050002

Accessories

- **Picture**
Connection parts
-



- **Connection parts:**

Please refer to our spare parts list for available tube adaptors.