

Now you can improve your Matrix™ Hood System!

Accutron's Personal Inhaler Plus (PIP+) Single Hood Sets incorporate a clear scavenging cone with a special thermoplastic elastomer construction to maintain good hood shape. The 100% latex free material minimises potential for allergic reactions in three comfortable sizes to accommodate all patient types.

The Benefits of the PIP+ Hood



By upgrading your existing Matrix™ Hood to the Accutron PIP+ System you benefit by adding a number of advantages to your existing tubing set.

- PIP+'s hood design promotes better access by being angled away from the oral area.

- The design of the PIP+ hoods creates both a comfortable and effective seal against the patient nasal region.


- Replacement of the hoods becomes easier with the PIP+ hood, with less effort required and minimal misalignment of the hood.

- The scavenging hood exhalation valve is included with every hood so the correct function is guaranteed

- PIP+'s clear scavenging cone allows the ability to monitor patients breathing, and verify the correct operation of the exhalation valve.

- The colours and assorted scents of the PIP+ Single Use Hoods combine to assist with overcoming patient anxiety, especially with children.

- Cost of the PIP+ hoods is competitive.



Original Matrix™ Hood Set
With Matrix scavenging cone

Matrix™ Set with PIP+ Upgrade
With adapter kit scavenging cone

PIP+ hood sets are available in Single Use Assorted Scented Types in Boxes of 24.

The coloured hoods offer a pleasant appearance with the scented flavours of the single use hoods that combine to assist with minimising any of the anxiety patient's experience with visiting the dentist.







Matrix™ ANS Tubing Set
with PIP+ Upgrade



All PIP+ hoods can easily be adapted to Matrix™ tubing sets, using an easy to fit adapter supplied free of charge*.



PIP+ Single Use Scavenging Hoods

	Code	Product	Features	List Price
	AC.33015-3	Large Variety Set	PIP+ Single Use Hood - Large Size Scented Variety Set - Mint - Peach - Vanilla Packet of 24	\$ 170+GST = \$ 187.00
	AC.33016-1	Medium Variety Set A	PIP+ Single Use Hood - Medium Size Scented Variety Set - Orange - Strawberry - Bubblegum Packet of 24	\$ 160+GST = \$ 176.00
	AC.33016-3	Medium Variety Set B	PIP+ Single Use Hood - Medium Size Scented Variety Set - Mint - Peach - Vanilla Packet of 24	\$ 160+GST = \$ 176.00
	AC.33017-1	Small Variety Set	PIP+ Single Use Hood - Small Size Scented Variety Set - Orange - Strawberry - Bubblegum Packet of 24	\$ 150+GST = \$ 165.00

PIP+ Autoclavable Scavenging Hoods

	Code	Product	Features	List Price
	AC.33018	Large Size Hood	Multi-Use Autoclavable Hood - Large Size Grey Colour One piece	\$ 170+GST = \$ 187.00
	AC.33019	Medium Size Hood	Multi-Use Autoclavable Hood - Medium Size Grey Colour One piece	\$ 160+GST = \$ 176.00
	AC.33020	Small Size Hood	Multi-Use Autoclavable Hood - Small Size Grey Colour One piece	\$ 150+GST = \$ 165.00

PIP+ Adapter Set for Matrix™ Tubing Circuits

	Code	Product	Features	List Price
	AC.35350	Adapter Set for Matrix Tubing Circuit	Adapter Set for PIP+ Single Hood to fit onto Matrix Tubing Circuit	No charge *
	AC.35300	Adapter Set for Matrix ANS Tubing Circuit	Adapter Set for PIP+ Single Hood to fit onto Matrix ANS Tubing Circuit - older version	No charge *

* No charge for Adapter Set when ordered with first order of PIP+ Single Use Hoods.